



DEVICE PARTS	DESCRIPTION	
	Device used for the filling of 100 ml crn® CRONO® Réservoir.	
	CLASSIFICATION	
	Class I medical device, according to Rule 12 of Annex IX of 93/42/EEC European directive and subsequent changes introduced by the 2007/47/EC European directive.	
	MANUFACTURER	
	CANÈ S.p.A. via Cuornè 42/a 10098 Rivoli (TO) - Italy Tel. +39-011-9574872 www.canespa.it mailbox@canespa.it	
	MARKING	
	ITALIAN DEVICE REPERTORY REGISTRATION	1772908
	ITALIAN DEVICE CLASSIFICATION (CND)	A0280
REF CANÈ	CRONO JET 100	
GTIN CODE	08050616619105	
PART N.	DESCRIPTION	MATERIAL
1	Reservoir attachment	ABS
2	O-ring	NBR
3	Rocker Switch	PA
4	Label	PET
5	Battery compartment	ABS


WARNING

- The **CRONOJET¹⁰⁰** was designed **esclusively** to fill 100ml CRN® **CRONO® Réservoir**.
- The use of the **CRONOJET¹⁰⁰** with devices different from 100ml CRN® **CRONO® Réservoir** can cause damage to both the syringe and the user.
- CANÈ S.p.A. disclaims all responsibility in the following cases:
 - the **CRONOJET¹⁰⁰** is used with other devices;
 - the **CRONOJET¹⁰⁰** is used with non-original reservoirs;
 - the **CRONOJET¹⁰⁰** is used out of its intended purpose.

PACKAGING	PCS	FEATURES	MATERIAL
Primary	1	Cardboard box 24 x 14 x 6 cm	White Cardboard

TECHNICAL SPECIFICATIONS

SPECIFIC	VALUE	NOTES
Supply voltage	$V_{al} = 3V$ dc	2 AA batteries
Intake pressure	$P_a = 0,4 \pm 0,1$ bar	Continuous mode
Outlet pressure	$P_s = 0,5 \pm 0,2$ bar	Continuous mode
Dimensions	124 x 71 x 56 mm	-
Weight	178 g	Battery included